



OBJECTIVE

Demonstrate the electric current generated by the motion of charged droplets of water

SUMMARY

Electric current arises due to an amount of charge being transported during a given interval of time. A flow of current can be simply illustrated with the help of charged droplets of water. In order to carry out the measurement, a burette and a Faraday cup connected to an electrometer will be used. The charge accumulated in the Faraday cup in a certain period of time is measured with the help of the voltage which drops across the capacitor. This allows the charge per droplet and the current to be determined.

EXPERIMENT PROCEDURE

- Measure the charge transferred to a Faraday cup by charged droplets of water dripping from a burette as a function of the time.
- Determine the current generated by the movement of the charged water droplets.
- Determine the charge on each droplet.

REQUIRED APPARATUS

Quantity	Description	Number
1	Electrometer (230 V, 50/60 Hz)	1001025 or
	Electrometer (115 V, 50/60 Hz)	1001024
1	Electrometer Accessories	1006813
1	Analogue Multimeter AM50	1003073
1	Burette, 10 ml	1018065
1	Constantan Wire 0.2 mm / 100 m	1000955
1	DC Power Supply 450 V (230 V, 50/60 Hz)	1008535 or
	DC Power Supply 450 V (115 V, 50/60 Hz)	1008534
1	Digital Multimeter P3340	1002785
1	Digital Stopwatch	1002811
1	Tripod Stand 150 mm	1002835
1	Stainless Steel Rod 1000 mm	1002936
2	Universal Clamp	1002830
1	Universal Jaw Clamp	1002833
1	Crocodile Clip 4 mm, Not Insulated	1002844
1	Set of 3 Safety Experiment Leads for Free Fall Apparatus	1002848
2	Pair of Safety Experimental Leads, 75cm, red/blue	1017718
1	Peleus ball, standard	1013392
1	Set of 10 Beakers, Low Form	1002872
Additionally recommended:		
1	3B NETlog™ (230 V, 50/60 Hz)	1000540 or
	3B NETlog™ (115 V, 50/60 Hz)	1000539
1	3B NETlab™	1000544



BASIC PRINCIPLES

Electric current arises due to an amount of charge being transported during a given interval of time. A flow of current can be simply illustrated with the help of charged droplets of water.

In this experiment a number of charged water droplets N drips at a constant rate of roughly one droplet per second from a burette into a Faraday cup connected to an electrometer and a capacitor. The charge Q accumulated in the Faraday cup causes the capacitor to charge up. The resulting voltage across the capacitor is observed and measured using an analogue multimeter for a certain period of time t . The high-resistance input of the operational amplifier in the electrometer ensures that the capacitor does not discharge via that path.

Observation of the analogue multimeter indicates that the voltage across the capacitor increases by about the same amount with every droplet collected in the Faraday cup, i.e. each of the droplets carries approximately the same charge:

$$(1) \quad q = \frac{Q}{N}$$

The current transported is given by

$$(2) \quad I = \frac{Q}{t}$$

As an option, the voltage across the capacitor can be recorded with the help of 3B NETlog™ and 3B NETlab™ as a function of time t and displayed in the form of a graph.

EVALUATION

The charge Q accumulated in the Faraday cup is determined by reading of the voltage U and calculating Q from that:

$$Q = C \cdot U \text{ where } C = 1 \text{ nF: capacitance of capacitor}$$

Using 3B NETlog™ and 3B NETlab™ the time characteristic $Q(t)$ can be measured. It is step-shaped whereby the individual steps mark the charge q that accumulates with each individual droplet per time interval Δt . The fact that each water droplet carries almost the same charge is reflected in the constant step height of the characteristic.




Fig. 1: Schematic illustrating the principle behind the measurement




Fig. 2: Accumulated charge Q as a function of time t