



NEW

## P22 SUTURE SKIN TRAINER


3B Scientific silicon suture pad perfectly replicates the anatomical structure of human skin layers including the epidermis, dermis and hypodermis to create a realistic haptic experience of touch for the user. Our objective is to assist students in acquiring the correct suturing technique through ongoing training.

The integration of a fabric net, both horizontally and vertically, within the silicone skin layers enhances durability and tear resistance, ensuring that the skin pad remains intact during suture training. Using the provided proprietary base, users can apply tension to the suture pad, resulting in a more realistic suturing experience. This tension replicates the alignment of the underlying collagen fibers. The skin pad allows users to make incisions of varying depth and location with a sharp scalpel. Each incision can be sutured repeatedly until a specific area on the pad becomes too perforated to support further suturing.




**Key Skills:**

- Subdermal suture
- Single interrupted suture
- Continuous suture
- Intracutaneous suture
- Suture stapling
- Z-Plasty suture



**Delivery content:**

- Suture pad (x2)
- Skin Pad Base



Skin Suture Set Light | 1024494



Skin Suture Set Dark | 1024495

**Replacement parts:**



Suture Pad Replacement  
10x light | 1024496



Suture Pad Replacement  
10x dark | 1024497

**Also available:**

Suture Practice Arm light  
1020904 | P101

Suture Practice Arm dark  
1023312 | P101D



MORE INFORMATION AVAILABLE ONLINE AT [3BSCIENTIFIC.COM](http://3BSCIENTIFIC.COM)



3B Scientific GmbH • Ludwig-Erhard-Straße 20 • 20459 Hamburg • Germany • [3bscientific.com](http://3bscientific.com)

Phone: +49 40 73966-0 • Fax: +49 40 73966-100 • E-mail: [info@3bscientific.com](mailto:info@3bscientific.com)